

## Experience Sharing Seminar for Food Testing Laboratories

# Testing of other harmful substances listed in the Harmful Substances in Food (Amendment) Regulation 2021 (Cap. 132AF)

8 December 2022

CHOI Lai-sheung, Timmy 蔡麗嫦  
Chemist / ACC Section  
Government Laboratory

Harmful substances in The Harmful Substances in Food (Amendment) Regulation 2021 (Cap. 132AF) mentioned in this presentation include

- Benzo[a]pyrene
- Erucic acid
- 3-MCPD
- Glycidyl fatty acid esters
- Partially hydrogenated oil
- (PHO)

Other information on testing methods are also available on CFS's website:

<https://www.cfs.gov.hk/harmfulsubstance/>

# Benzo[a]pyrene

## Oil and fat

- Weigh sample
- Add  $d_{12}$ -BaP and cyclohexane
- Proceed to MIP-SPE cleanup

## Infant formula

- Weigh sample
- Add  $d_{12}$ -BaP, EtOH, water and 25%  $\text{NH}_3$

## Fat extraction

- Warm at 65 °C.
- Cool to r.t. Add EtOH and transfer to Mojonnier flask
- Add **Et<sub>2</sub>O** and **petroleum ether**, shake
- Centrifuge and collect the upper layer
- Repeat extraction
- Combine the organic layers and evaporate solvent
- Reconstitute with cyclohexane
- Proceed to MIP-SPE cleanup

## Molecularly-imprinted polymer (MIP)-SPE cleanup

- Condition **MIP-SPE** with cyclohexane
- Load sample onto MIP-SPE column
- Wash with cyclohexane
- Elute with EA
- Evaporate eluate and reconstitute with pentane

## Further cleanup

- Load sample onto **deactivated silica gel column**
- Collect eluate
- Elute with pentane
- Evaporate eluate to about 0.2 mL
- Add  $d_{12}$ -perylene as injection standard

## GC-MS/MS analysis

Scope:

*Oil, fat, infant and follow-up formula*

For reference:

*GB5009.265-2021*

# Benzo[a]pyrene

## GC-MS/MS analysis

### GC settings

Column	:	DB-35ms (30 m × 0.25 mm, 0.25 μm)
Injector temperature	:	280 °C
Injector mode	:	Splitless
Injection volume	:	1 μL
Column flow	:	Helium at 1.2 mL/min
Oven programme	:	150 °C (hold for 1 min), 20 °C/min to 320 °C (hold for 15 min)

### MS settings

MS system	:	TSQ 8000 Evo
Ionization mode	:	EI
Transfer line temperature	:	280 °C
Ion source temperature	:	300 °C

Analytes	Retention time (min)	Precursor ion (m/z)	Product ion (m/z)	Collision energy (eV)
d <sub>12</sub> -BaP	12.9	264	260	25
BaP	13.0	252	250	30
			251	15
d <sub>12</sub> -Perylene	13.2	264	260	25

# Erucic acid

## *Scope:*

*Oil, fat and food to which oil and/or fat has been added*

## *For reference:*

*AOAC 922.06 (fat extraction);*

*ISO 12966-2:2017, AOAC 963.22, AOAC 969.33*

## Oil and fat

- Weigh sample into RB flask

## Transmethylation

- Add **methanolic NaOH**
- Reflux until fat droplets disappear
- Add **12.5% (w/v) BF<sub>3</sub> in MeOH**, boil
- Add heptane and boil
- Remove from heat
- Add sat. NaCl until the heptane layer float to the neck of flask
- Transfer out the heptane layer and dry with anhy. Na<sub>2</sub>SO<sub>4</sub>

## GC-FID analysis

# Erucic acid

## GC-FID analysis

Column	:	HP-DB-FFAP (30 m × 0.53 mm, 1 μm)
Injector temperature	:	250 °C
Injection mode	:	Splitless
Injection volume	:	1 μL
Column flow	:	Helium at 7 mL/min
Column temperature	:	170 °C (hold for 1 min), 10 °C/min to 210 °C (hold for 6 min), 10 °C/min to 240 °C (hold for 9 min).
Detector temperature	:	250 °C
Make up gas and flow	:	Helium, 15 mL/min
H <sub>2</sub> flow	:	40 mL/min
Air flow	:	400 mL/min

*For reference:  
ISO 12966-4:2015*

# 3-MCPD

*Scope:  
Solid and liquid condiments containing  
acid hydrolysed vegetable proteins*

*For reference:  
BS EN ISO 14573:2004*

## Condiments

- Add  $d_5$ -3-MCPD into sample

## Removal of fat

- Add **heptane** and 5 M NaCl
- Mix and separate the layers
- Discard the heptane layer

## Supported Liquid Extraction (SLE) Cleanup

- Mix aqueous layer with **diatomaceous earth-based solid phase**
- Load into column packed with the same solid phase
- Elute with **Et<sub>2</sub>O**
- Evaporate to ~2 mL, dry with an. Na<sub>2</sub>SO<sub>4</sub>
- Transfer out the dried eluate, add 0.1 mL isooctane
- Evaporate solvent to about 0.1 mL

## Derivatisation

- Add isooctane and heptafluorobutyric imidazole (**HFBI**)
- Incubate at 70 °C
- Add water, mix
- Separate the isooctane layer, dry with anhy. Na<sub>2</sub>SO<sub>4</sub>
- Collect the isooctane layer

## GC-MS analysis

# 3-MCPD

## GC-MS analysis

### GC settings

Column	:	HP-5MS (30 m × 0.25 mm, 0.25 μm)
Injector temperature	:	250 °C
Injector mode	:	Splitless
Injection volume	:	1 μL
Column flow	:	Helium at 1 mL/min (constant flow)
Oven programme	:	58 °C (hold for 3 min), 1 °C/min to 65 °C, 2 °C/min to 77 °C, 50 °C/min to 300 °C (hold for 10 min).



# 3-MCPD

## GC-MS analysis

### MS settings

MS system	:	Shimadzu QP2020 NX
Ionization mode	:	EI
Transfer line temperature	:	280 °C
Source temperature	:	230 °C

MS system	:	Q Exactive GC Orbitrap
Ionization mode	:	EI (70 eV)
Transfer line temperature	:	280 °C
Ion source temperature	:	280 °C
Resolution	:	120,000
AGC target	:	5e6
Maximum injection time	:	Auto
Microscan	:	1
Scan range	:	m/z 50 - 500

Analytes (derivatised product)	Retention time (min)	Ions monitored		
		Quantitation ion	Qualifier ions	
d <sub>5</sub> -3-MCPD (+ 2HFB)	14.4	257.03451	-	-
3-MCPD (+ 2HFB)	14.6	274.97067	452.98023	288.98616
<b>2-MCPD</b> (+ 2HFB)	14.9	74.99960	76.99665	288.98616

# Glycidyl fatty acid esters

Glycidyl fatty acid esters expressed as glycidol

## Solid and liquid infant formula

- Add  $d_5$ -glycidyl ester

## Fat extraction

- Add **EA** and water
- Sonicate at 50 °C
- Separate and collect the EA layer
- Repeat EA extraction
- Combine EA layers. Evaporate to dryness
- Reconstitute with THF

## Conversion of GE to 3-MBPDE

- Add **acidic NaBr**
- Incubate at 50 °C
- Add 0.6% (w/v)  $\text{Na}_2\text{CO}_3$
- Extract with heptane
- Evaporate the heptane layer
- Reconstitute with THF

## Methanolysis of 3-MBPDE to 3-MBPD

- Add **1.8% (v/v)  $\text{H}_2\text{SO}_4$  in MeOH**
- Incubate at 40 °C
- Add sat.  $\text{Na}_2\text{CO}_3$ , vortex.
- Evaporate solvent to 1 mL
- Add 20% (w/v)  $\text{Na}_2\text{SO}_4$  and heptane
- Vortex. Remove the organic layer
- Extract the aqueous layer with EA twice

## Derivatisation

- Add 0.1 mL isooctane to the combined EA layer
- Evaporate to 0.1 mL
- Add isooctane, **HFBI** and anhy.  $\text{Na}_2\text{SO}_4$
- Incubate at 70 °C
- Add water, vortex
- Collect the isooctane layer

## GC-MS analysis

*Scope:*

*Infant formula and follow-up formula*

*For reference:*

*AOAC 2018.03*

# Glycidyl fatty acid esters

Glycidyl fatty acid esters expressed as glycidol

GC-MS analysis

## GC settings

Same as the analysis of 3-MCPD.

## MS settings

MS system	:	Q Exactive GC Orbitrap
Ionization mode	:	EI (70 eV)
Transfer line temperature	:	280 °C
Ion source temperature	:	280 °C
Resolution	:	120,000
AGC target	:	5e6
Maximum injection time	:	Auto
Microscan	:	1
Scan range	:	m/z 50 - 500

Analytes (derivatised product)	Retention time (min)	Ions monitored		
		Quantitation ion	Qualifier ions	
d <sub>5</sub> -3-MBPD (+ 2HFB)	18.6	257.03451	-	-
3-MBPD (+ 2HFB)	18.7	120.94747	118.94925	253.00940

# Melamine

## Without SPE cleanup

- Weigh sample
- Make up with 49:49:2 (v/v/v) MeCN / water / 1 M HCl
- Transfer out 5 mL
- Wash with CH<sub>2</sub>Cl<sub>2</sub>
- Collect the upper aqueous layer
- Back wash the CH<sub>2</sub>Cl<sub>2</sub> layer with 0.1 M HCl
- Collect and combine the aqueous layers
- Make up with water
- Dilute with MeCN, filter
- Add <sup>13</sup>C<sub>3</sub><sup>15</sup>N<sub>3</sub>-melamine

## LC-MS/MS analysis

## With SPE cleanup

- Add <sup>13</sup>C<sub>3</sub><sup>15</sup>N<sub>3</sub>-melamine into sample
- Make up with 49:49:2 (v/v/v) MeCN / water / 1 M HCl
- Transfer out 5 mL
- Wash with CH<sub>2</sub>Cl<sub>2</sub>
- Collect the upper aqueous layer
- Back wash the CH<sub>2</sub>Cl<sub>2</sub> layer with 0.1 M HCl
- Collect and combine the aqueous layers
- Make up with water

## SPE cleanup

- Condition **MCX column** with MeOH and water
- Load aqueous sample
- Wash with 0.1 M HCl and MeOH
- Dry the column
- Elute with 5% ammonia in MeOH
- Remove solvent and reconstitute with 95:5 (v/v).0.1% formic acid in MeCN / 20 mM NH<sub>4</sub>OAc

## LC-MS/MS analysis

### Scope:

*Infant formula and follow-up formula, milk and food*

For reference:

[https://www.govtlab.gov.hk/doc/our\\_work/texchange/melamine.pdf](https://www.govtlab.gov.hk/doc/our_work/texchange/melamine.pdf)

# Melamine

## LC-MS/MS analysis

### LC settings

Column	:	Waters Altantis HILIC column (100 mm x 2.1 mm, 3 $\mu$ m particle size)
Mobile phase	:	A: 20 mM Ammonium acetate in water B: Acetonitrile
Flow rate	:	0.5 mL/min
Injection volume	:	5 $\mu$ L
Solvent programme	:	

Time (min)	Solvent A (% v/v)	Solvent B (% v/v)
0.0	5	95
3.8	95	5
6.0	95	5
6.5	5	95
13.5	5	95

# Melamine

## LC-MS/MS analysis

### MS settings

MS/MS system	:	Sciex Qtrap 5500
Ionization mode	:	+ ESI
Temperature (TEM)	:	550 °C
Curtain Gas (CUR)	:	25 psi
Collision Gas (CAD)	:	Medium
Ion Spray Voltage (IS)	:	5500 V
Ion Source Gas 1 (GS1)	:	55 psi
Ion Source Gas 2 (GS2)	:	45 psi

Analyte	Retention time (min)	Precursor ion (m/z)	Product ion (m/z)	Declustering Potential (DP)	Entrance Potential (EP)	Collision Energy (CE)	Collision Energy (CXP)
Melamine	2.94	127	85	81 V	10 V	25 V	22 V
			68	81 V	10 V	39 V	30 V
			60	81 V	10 V	27 V	14 V
Melamine- <sup>13</sup> C <sub>3</sub> <sup>15</sup> N <sub>3</sub>	2.94	133	89	96 V	10 V	25 V	16 V

# Partially hydrogenated oil (PHO)

## Cap132AF Amendment Regulation

Partially hydrogenated oil (PHO) means any oil or fat that has undergone the process of hydrogenation but is not fully saturated as a result of that process.

部分氫化油指有經過氫化的過程但並無因為該過程而完全飽和的油或脂肪。

- PHO is the major source of industrially produced trans fatty acids (IP-TFA).  
部分氫化油是工業生產的反式脂肪酸的主要來源。
- Natural sources of dietary trans fatty acids are available from ruminal source and refined oils.  
天然的膳食反式脂肪酸可於反芻動物和精煉油中獲得。
- Need to differentiate whether the detected TFA is from industrial source.  
需要區分檢測到的反式脂肪酸是否來自工業來源。

# IP-TFA

## EU's analytical approach

(Commission Regulation (EU) 2019/649)

1. **Total fat**
2. **Butyric acid** ( $C_{4:0}$ ) (g/100 g fat)
3. **Total TFA** (g/100 g fat): Sum of four non-conjugated TFAs –
  - hexadecenoic acid (trans- $C_{16:1}$ )
  - octadecenoic acid (trans- $C_{18:1}$ )
  - octadecadienoic acid (trans- $C_{18:2}$ )
  - octadecatrienoic acid (trans- $C_{18:3}$ )
4. **Conjugated linoleic acid (CLA)** (9-cis,11-trans- $C_{18:2}$ )

European Commission

JRC TECHNICAL REPORT

Analytical approach for checking the compliance of fats and oils against the regulated limit for industrial trans fatty acids (Commission Regulation (EU) 2019/649)


2021

Oleic acid (cis-D9)

Elaidic acid (trans-D9)

Joint Research Centre

EUR 30767 EN





## Analysis of total fat

### General foods

- Add EtOH to sample. Mix well

### Acid hydrolysis

- Add 8.3 M **HCl**. Mix well
- Shake in water bath at 70–80 °C
- Add EtOH and cool to r.t.

### Extraction of fat

- Transfer content to Mojonnier flask
- Add **Et<sub>2</sub>O** and shake
- Add **petroleum ether** and shake
- Centrifuge
- Collect the ether layer
- Repeat the extraction with Et<sub>2</sub>O and petroleum ether twice
- Combine the ether layer and remove solvent
- Dry in oven at 100 °C. Stand in desiccator

*For reference:*  
AOAC 922.06

### Weighing

# IP-TFA

## Analysis of fatty acids

- Weigh sample
- Add pyrogalllic acid and triundecanoin (C11:0 triglyceride)
- Add EtOH. Mix well.

### For dairy products

- Add H<sub>2</sub>O. Mix.
- Add **NH<sub>4</sub>OH** (58% w/w). Mix.
- Shake in water bath at 70–80 °C
- Keep solution basic to phenolphthalein indicator using NH<sub>4</sub>OH.
- Proceed to fat extraction.

### For cheese and formula products

- Add H<sub>2</sub>O. Mix.
- Add **NH<sub>4</sub>OH**. Mix.
- Shake in water bath at 70–80 °C.
- Add 12 M HCl and keep in boiling water bath.
- Proceed to fat extraction.

### For other foods

- Add 8.3 M **HCl**. Mix.
- Shake in water bath at 70–80 °C.
- Proceed to fat extraction.

## Extraction of fat

- Remove sample from water bath. Cool to r.t.
- Transfer content to Mojonnier flask
- Add **Et<sub>2</sub>O** and shake
- Add **petroleum ether** and shake
- Centrifuge. Collect the ether layer.
- Repeat the extraction with Et<sub>2</sub>O and petroleum ether twice.
- Filter and dry with Na<sub>2</sub>SO<sub>4</sub>.
- Evaporate solvent to about 2 mL

## Methylation

- Add **7% BF<sub>3</sub> in methanol** and toluene to sample.
- Heat at 100 °C
- Cool to r.t. and add H<sub>2</sub>O, hexane and Na<sub>2</sub>SO<sub>4</sub>. Shake.
- Collect the top layer and dry with Na<sub>2</sub>SO<sub>4</sub>.

## GC-FID analysis

*For reference:*  
AOAC 996.06

# IP-TFA

## GC-FID analysis

Column	:	SP-2560 (100 m × 0.25 mm, 0.2 μm)
Injector temperature	:	260 °C
Injection mode	:	Split (50:1)
Injection volume	:	1 μL
Column flow	:	Helium at 0.7 mL/min
Column temperature	:	100 °C (hold for 4 min), 3 °C/min to 240 °C (hold for 20 min).
Detector temperature	:	250 °C
Make up gas and flow	:	Helium, 45 mL/min
H <sub>2</sub> flow	:	40 mL/min
Air flow	:	400 mL/min

# Reference materials

## Certified Reference Materials

Benzo[a]pyrene	– NMIJ
Erucic acid	– Sigma-Aldrich (as mixed standard)
Melamine	– Supelco

## Reference Materials

Erucic acid	– Santa Cruz Biotechnology, Sigma-Aldrich
3-MCPD	– Sigma-Aldrich, TCI
2-MCPD	– Santa Cruz Biotechnology, TRC
Glycidyl palmitate	– TRC, Wako
Fatty acids	– Nu-Chek-Prep, Sigma-Aldrich

## Isotope Labelled Internal Standards

$d_{12}$ -BaP	– Cambridge Isotope Laboratory
$d_{12}$ -Perylene	– Cambridge Isotope Laboratory
$d_5$ -3-MCPD	– Cambridge Isotope Laboratory
$d_5$ -glycidyl palmitate	– Medical Isotopes
$^{13}\text{C}_3$ $^{15}\text{N}_3$ -melamine	– Cambridge Isotope Laboratory
Triundecanoin	– Nu-Chek-Prep

The End  
Thank you